

CSDE1 Rabbit mAb

Catalog No: #60252

Package Size: #60252-1 50ul #60252-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CSDE1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	CSDE1 Antibody detects endogenous levels of total CSDE1
Immunogen Description	A synthesized peptide derived from human CSDE1
Other Names	Csde1; NRU; UNR;
Accession No.	Uniprot:O75534
Calculated MW	89kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Background

RNA-binding protein. Required for internal initiation of translation of human rhinovirus RNA. May be involved in translationally coupled mRNA turnover.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.