

## TMF Rabbit mAb

Catalog No: #60258

Package Size: #60258-1 50ul #60258-2 100ul


## Description

Product Name	TMF Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human
Specificity	TMF Antibody detects endogenous levels of total TMF
Immunogen Description	A synthesized peptide derived from human TMF
Other Names	ARA160; TMF; Tmf1;
Accession No.	Uniprot:P82094
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200

## Images



Western blot analysis of TMF expression in HeLa cell lysate.

## Background

Potential coactivator of the androgen receptor. Mediates STAT3 degradation. May play critical roles in two RAB6-dependent retrograde transport processes: one from endosomes to the Golgi and the other from the Golgi to the ER. This protein binds the HIV-1 TATA element and inhibits transcriptional activation by the TATA-binding protein (TBP).

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.