

Recombinant Human Protein FAM110B(FAM110B)

Catalog No: #AP70290

Package Size: #AP70290-1 20ug #AP70290-2 100ug #AP70290-3 1mg

Description

Product Name	Recombinant Human Protein FAM110B(FAM110B)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-370aaSequence Info:Full Length
Accession No.	Q8TC76
Calculated MW	44.7 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	MPTETLQTGSMVKPVSPAGTFTSAVPLRILNKGPDYFRRQAEPNPKRLSAVERLEADKAKYVKSQEVINAKQE PVKPAVLAKPPVCPAAKRALGSPTLKVFGNHAKTESGVQRENKLEILKNIINSSEGSSSGSGHKHSSRNWPP HRSEATDLHRHSFAESLKVYPTQGRRSPQEGGSHVGRRLLEQSAESFLHVSHSSDIRKVTSVKPLKAIPCSS SAPPLPPKPKIAAIASMKSPADPVEPACGVSRRLQSRKSDLSDRYFRVDADVERFFNYCGLDPEELENLG MENFARANSDIISLNFASMISSDCEQSQSDNSDLRNDDSANRVPYGISAIERNARIKWLYSIKQARESQKV SHV
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May be involved in tumor progression.

References

The consensus coding sequences of human breast and colorectal cancers.Sjoebloom T., Jones S., Wood L.D., Parsons D.W., Lin J., Barber T.D., Mandelker D., Leary R.J., Ptak J., Silliman N., Szabo S., Buckhaults P., Farrell C., Meeh P., Markowitz S.D., Willis J., Dawson D., Willson J.K.V. , Gazdar A.F., Hartigan J., Wu L., Liu C., Parmigiani G., Park B.H., Bachman K.E., Papadopoulos N., Vogelstein B., Kinzler K.W., Velculescu V.E.Science 314:268-274(2006)Research Topic:Cancer

Note: This product is for in vitro research use only and is not intended for use in humans or animals.