

Recombinant Pig Interleukin-33(IL33),partial

Catalog No: #AP70441

Package Size: #AP70441-1 20ug #AP70441-2 100ug #AP70441-3 1mg

Description

Product Name	Recombinant Pig Interleukin-33(IL33),partial
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:123-276aaSequence Info:Partial
Accession No.	M5B263
Calculated MW	44.6 kDa
Tag Info	N-terminal GST-tagged
Target Sequence	EHSASLSTYNDQYITFAFEDGSYEIYVEDLRKDQEKDKVLLRYYDSQIPSETDGGGDHRKLMVNLSP TKDKD FLLHANSKEHSVELQKCENPLPEQAFFVLHEQPSKCVSFECKSHPGVFLGVKNNQLALIKLGEHPEDSNRENT TFKLSNLM
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

"cDNA cloning and expression of porcine IL-33."Shimazu T., Tohno M., Kawai Y., Saito T., Kitazawa H.Submitted (FEB-2007) Research
Topic:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.