

Recombinant Human Sorting nexin-24(SNX24)

Catalog No: #AP70794

Package Size: #AP70794-1 20ug #AP70794-2 100ug #AP70794-3 1mg

Description

Product Name	Recombinant Human Sorting nexin-24(SNX24)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-169aaSequence Info:Full Length
Accession No.	Q9Y343
Calculated MW	35.8 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	MEVYIPSFYEEEDLGRGYTVFKIEVLMNGRKHFVEKRYSEFHALHKKLKKCIKTPEIPSKHVRNWVPKVLEQR RQGLETYLQAVILENEELPKLFLDFLNVRHLPSPKAESCGSFDETESEESSKLSHQPVLLFLRDPYVLPAAASD FPNVVIEGVLHGIFYPHLQPR
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May be involved in several stages of intracellular trafficking.

References

N-terminal acetylome analyses and functional insights of the N-terminal acetyltransferase NatB.Van Damme P., Lasa M., Polevoda B., Gazquez C., Elosegui-Artola A., Kim D.S., De Juan-Pardo E., Demeyer K., Hole K., Larrea E., Timmerman E., Prieto J., Arnesen T., Sherman F., Gevaert K., Aldabe R.Proc. Natl. Acad. Sci. U.S.A. 109:12449-12454(2012) Research Topic:Transport

Note: This product is for in vitro research use only and is not intended for use in humans or animals.