

Recombinant Rat Sortilin(Sort1),partial

Catalog No: #AP70798

Package Size: #AP70798-1 20ug #AP70798-2 100ug #AP70798-3 1mg

Description

Product Name	Recombinant Rat Sortilin(Sort1),partial
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:610-754aaSequence Info:Partial
Other Names	Glycoprotein 110 ;Gp110Neurotensin receptor 3 ;NTR3
Accession No.	O54861
Calculated MW	32.5 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	CEENDYTTWLAHSTDPGDYKDGCGILGYKEQFLRLRKSSVCQNGRDYVVAQKPSICPCSLEDFLCDFGYFRPE NASECVEQPELKGHELEFCLYGKEEHLTTNGYRKIPGDRCQGGMNPAREVKDLKKKCTSNFLNPKKQNSKSS S
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Functions as a sorting receptor in the Golgi compartment and as a clearance receptor on the cell surface. Required for protein transport from the Golgi apparatus to the lysosomes by a pathway that is independent of the mannose-6-phosphate receptor (M6PR). Also required for protein transport from the Golgi apparatus to the endosomes. Promotes neuronal apoptosis by mediating endocytosis of the proapoptotic precursor forms of BDNF (proBDNF) and NGFB (proNGFB). Also acts as a receptor for neurotensin. May promote mineralization of the Extracellular domain matrix during osteogenic differentiation by scavenging Extracellular domain LPL. Probably required in adipocytes for the formation of specialized storage vesicles containing the glucose transporter SLC2A4, GLUT4 (GLUT4 storage vesicles, or GSVs). These vesicles provide a stable pool of SLC2A4 and confer increased responsiveness to insulin. May also mediate transport from the endoplasmic reticulum to the Golgi.

References

Sortilin controls intracellular sorting of brain-derived neurotrophic factor to the regulated secretory pathway.Chen Z.-Y., Ieraci A., Teng H., Dall H., Meng C.-X., Herrera D.G., Nykjaer A., Hempstead B.L., Lee F.S.J. Neurosci. 25:6156-6166(2005) Research Topic:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.