

Recombinant Human AQP4/Aquaporin-4 Protein, N-His-KSI

Catalog No: #AP97687

Package Size: #AP97687-1 50ug #AP97687-2 100ug #AP97687-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AQP4/Aquaporin-4 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Aquaporin-4, AQP-4, Mercurial-insensitive water channel, MIWC, WCH4, AQP4
Accession No.	P55087
Target Species	Homo sapiens (Human)
Calculated MW	25.71 kDa
Target Length	[Thr137-Thr155] & [Cys253-Val323]
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.