

Recombinant Human DST/Dystonin Protein, N-His

Catalog No: #AP97695

Package Size: #AP97695-1 50ug #AP97695-2 100ug #AP97695-3 1mg

Description

Product Name	Recombinant Human DST/Dystonin Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Dystonin, 230 kDa bullous pemphigoid antigen, 230/240 kDa bullous pemphigoid antigen, Bullous pemphigoid antigen 1, BPA, Bullous pemphigoid antigen, Dystonia musculorum protein, Hemidesmosomal plaque protein, DST, BP230, BP240, BPAG1, DMH, DT, KIAA0728
Accession No.	Q03001
Target Species	Homo sapiens (Human)
Calculated MW	28.97 kDa
Target Length	Lys580-Glu803
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.