

Recombinant Human PHF5A Protein, N-His-SUMO & C-Strep

Catalog No: #AP97752

Package Size: #AP97752-1 50ug #AP97752-2 100ug #AP97752-3 1mg

Description

Product Name	Recombinant Human PHF5A Protein, N-His-SUMO & C-Strep
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	PHD finger-like domain-containing protein 5A, PHD finger-like domain protein 5A, Splicing factor 3B-associated 14 kDa protein, SF3b14b, PHF5A
Accession No.	Q7RTV0
Target Species	Homo sapiens (Human)
Calculated MW	23.85 kDa
Target Length	Arg21-Arg110
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.