

Recombinant Human WWP1 Protein, N-His

Catalog No: #AP97759

Package Size: #AP97759-1 50ug #AP97759-2 100ug #AP97759-3 1mg

Description

Product Name	Recombinant Human WWP1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	NEDD4-like E3 ubiquitin-protein ligase WWP1, 2.3.2.26, Atrophin-1-interacting protein 5, AIP5, HECT-type E3 ubiquitin transferase WWP1, TGIF-interacting ubiquitin ligase 1, Tiul1, WW domain-containing protein 1, WWP1
Accession No.	Q9H0M0
Target Species	Homo sapiens (Human)
Calculated MW	16.28 kDa
Target Length	Arg17-Gly141
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.