

Recombinant Human CHCHD6 Protein, N-His

Catalog No: #AP97760

Package Size: #AP97760-1 50ug #AP97760-2 100ug #AP97760-3 1mg

Description

Product Name	Recombinant Human CHCHD6 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	MICOS complex subunit MIC25, Coiled-coil-helix cristae morphology protein 1, Coiled-coil-helix-coiled-coil-helix domain-containing protein 6, CHCHD6, CHCM1, MIC25
Accession No.	Q9BRQ6
Target Species	Homo sapiens (Human)
Calculated MW	20.43 kDa
Target Length	Gln79-Gly235
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.