

Recombinant Human FOXO6 Protein, N-His-SUMO

Catalog No: #AP97761

Package Size: #AP97761-1 50ug #AP97761-2 100ug #AP97761-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FOXO6 Protein, N-His-SUMO
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Forkhead box protein O6, FOXO6
Accession No.	A8MYZ6
Target Species	Homo sapiens (Human)
Calculated MW	21.65 kDa
Target Length	Ser92-Gly172
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.