

Recombinant Human TSPYL2 Protein, N-His

Catalog No: #AP97762

Package Size: #AP97762-1 50ug #AP97762-2 100ug #AP97762-3 1mg

Description

Product Name	Recombinant Human TSPYL2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	NP79, Nuclear protein of 79 kDa, CTCL-associated antigen se20-4, Cutaneous T-cell lymphoma-associated antigen se20-4, TSPY-like protein 2, DENTT, Differentially-expressed nucleolar TGF-beta1 target protein, Testis-specific Y-encoded-like protein 2, CDA1, Cell division autoantigen 1, TSPX, TSPYL2
Accession No.	Q9H2G4
Target Species	Homo sapiens (Human)
Calculated MW	21.28 kDa
Target Length	Pro263-Arg418
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.