

Recombinant Human PPARGC1B Protein, N-His

Catalog No: #AP97769

Package Size: #AP97769-1 50ug #AP97769-2 100ug #AP97769-3 1mg

Description

Product Name	Recombinant Human PPARGC1B Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Peroxisome proliferator-activated receptor gamma coactivator 1-beta, PGC-1-beta, PPAR-gamma coactivator 1-beta, PPARGC-1-beta, PGC-1-related estrogen receptor alpha coactivator, PPARGC1B, PERC, PGC1, PGC1B, PPARGC1
Accession No.	Q86YN6
Target Species	Homo sapiens (Human)
Calculated MW	18.32 kDa
Target Length	Ser886-His1023
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.