

Recombinant Human OOEP Protein, N-His

Catalog No: #AP97783

Package Size: #AP97783-1 50ug #AP97783-2 100ug #AP97783-3 1mg

Description

Product Name	Recombinant Human OOEP Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Oocyte-expressed protein homolog, KH homology domain-containing protein 2, Oocyte- and embryo-specific protein 19, hOEP19, OOEP, C6orf156, KHDC2, OEP19
Accession No.	A6NGQ2
Target Species	Homo sapiens (Human)
Calculated MW	17.18 kDa
Target Length	Gln23-Ala149
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.