

Recombinant Human FANCM Protein, N-His

Catalog No: #AP97792

Package Size: #AP97792-1 50ug #AP97792-2 100ug #AP97792-3 1mg

Description

Product Name	Recombinant Human FANCM Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Fanconi anemia group M protein, Protein FACM, 3.6.4.13, ATP-dependent RNA helicase FANCM, Fanconi anemia-associated polypeptide of 250 kDa, FAAP250, Protein Hef ortholog, FANCM, KIAA1596
Accession No.	Q8IYD8
Target Species	Homo sapiens (Human)
Calculated MW	41.44 kDa
Target Length	Ser281-His620
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.