

Recombinant Human IKZF3 Protein, N-His

Catalog No: #AP97801

Package Size: #AP97801-1 50ug #AP97801-2 100ug #AP97801-3 1mg

Description

Product Name	Recombinant Human IKZF3 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Zinc finger protein Aiolos, Ikaros family zinc finger protein 3, IKZF3, ZNFN1A3
Accession No.	Q9UKT9
Target Species	Homo sapiens (Human)
Calculated MW	30.26 kDa
Target Length	Glu84-Thr326
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.