

Recombinant Mouse CD45/PTPRC Protein, N-His

Catalog No: #AP97805

Package Size: #AP97805-1 50ug #AP97805-2 100ug #AP97805-3 1mg

Description

Product Name	Recombinant Mouse CD45/PTPRC Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Receptor-type tyrosine-protein phosphatase C, 3.1.3.48, Leukocyte common antigen, L-CA, Lymphocyte antigen 5, Ly-5, T200, CD45, Ptpcr, Ly-5
Accession No.	P06800
Target Species	Mus musculus (Mouse)
Calculated MW	43.04 kDa
Target Length	Thr202-Phe563
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.