

Recombinant Mouse CD20/MS4A1 Protein, N-GST & C-His

Catalog No: #AP97808

Package Size: #AP97808-1 50ug #AP97808-2 100ug #AP97808-3 1mg

Description

Product Name	Recombinant Mouse CD20/MS4A1 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	B-lymphocyte antigen CD20, B-cell differentiation antigen Ly-44, Lymphocyte antigen 44, Membrane-spanning 4-domains subfamily A member 1, CD20, Ms4a1, Cd20, Ly-44, Ms4a2
Accession No.	P19437
Target Species	Mus musculus (Mouse)
Calculated MW	33.98 kDa
Target Length	Leu133-Ser182
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.