

Recombinant Mouse CD172a/SIRPA Protein, N-His

Catalog No: #AP97818

Package Size: #AP97818-1 50ug #AP97818-2 100ug #AP97818-3 1mg

Description

Product Name	Recombinant Mouse CD172a/SIRPA Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Tyrosine-protein phosphatase non-receptor type substrate 1, SHP substrate 1, SHPS-1, Brain Ig-like molecule with tyrosine-based activation motifs, Bit, CD172 antigen-like family member A, Inhibitory receptor SHPS-1, MyD-1 antigen, Signal-regulatory protein alpha-1, Sirp-alpha-1, mSIRP-alpha1, p84, CD172a, Sirpa, Bit, Myd1, Ptpns1, Shps1, Sirp
Accession No.	P97797
Target Species	Mus musculus (Mouse)
Calculated MW	37.66 kDa
Target Length	Gly31-Gly350
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.