

Recombinant Mouse CD279/PDCD1/PD1 Protein, N-His

Catalog No: #AP97822

Package Size: #AP97822-1 50ug #AP97822-2 100ug #AP97822-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse CD279/PDCD1/PD1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Programmed cell death protein 1, Protein PD-1, mPD-1, CD279, Pdcd1, Pd1
Accession No.	Q02242
Target Species	Mus musculus (Mouse)
Calculated MW	15.46 kDa
Target Length	Gly30-Glu146
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.