

Recombinant Mouse VWF Protein, N-GST & C-His

Catalog No: #AP97854

Package Size: #AP97854-1 50ug #AP97854-2 100ug #AP97854-3 1mg

Description

Product Name	Recombinant Mouse VWF Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	von Willebrand factor, vWF, von Willebrand antigen 2, von Willebrand antigen II, Vwf
Accession No.	Q8CIZ8
Target Species	Mus musculus (Mouse)
Calculated MW	36.30 kDa
Target Length	Asp1596-Thr1668
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.