

Recombinant Mouse ISG15 Protein, N-His

Catalog No: #AP97857

Package Size: #AP97857-1 50ug #AP97857-2 100ug #AP97857-3 1mg

Description

Product Name	Recombinant Mouse ISG15 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Ubiquitin-like protein ISG15, Interferon-induced 15 kDa protein, Interferon-induced 17 kDa protein, IP17, Ubiquitin cross-reactive protein, Isg15, G1p2, Ucrp
Accession No.	Q64339
Target Species	Mus musculus (Mouse)
Calculated MW	19.15 kDa
Target Length	Met1-Arg151
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.