

## Recombinant Mouse RANBP2 Protein, N-His

Catalog No: #AP97893

Package Size: #AP97893-1 50ug #AP97893-2 100ug #AP97893-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                |                                                                                                                                                                           |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant Mouse RANBP2 Protein, N-His                                                                                                                                   |
| Host Species   | E. coli                                                                                                                                                                   |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                           |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                           |
| Other Names    | E3 SUMO-protein ligase RanBP2, 2.3.2.-, Ran-binding protein 2, RanBP2, Ranbp2                                                                                             |
| Accession No.  | Q9ERU9                                                                                                                                                                    |
| Target Species | Mus musculus (Mouse)                                                                                                                                                      |
| Calculated MW  | 21.06 kDa                                                                                                                                                                 |
| Target Length  | Tyr1840-Ala2000                                                                                                                                                           |
| Tag Info       | His-tag                                                                                                                                                                   |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.