

Recombinant Mouse HSD11B2 Protein, N-His

Catalog No: #AP97901

Package Size: #AP97901-1 50ug #AP97901-2 100ug #AP97901-3 1mg

Description

Product Name	Recombinant Mouse HSD11B2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	11-beta-hydroxysteroid dehydrogenase type 2, 11-DH2, 11-beta-HSD2, 1.1.1.-, Corticosteroid 11-beta-dehydrogenase isozyme 2, NAD-dependent 11-beta-hydroxysteroid dehydrogenase, Hsd11b2, Hsd11k
Accession No.	P51661
Target Species	Mus musculus (Mouse)
Calculated MW	29.35 kDa
Target Length	Thr82-Ala331
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.