

Recombinant Mouse MAPK7 Protein, N-His

Catalog No: #AP97907

Package Size: #AP97907-1 50ug #AP97907-2 100ug #AP97907-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse MAPK7 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Mitogen-activated protein kinase 7, MAP kinase 7, MAPK 7, 2.7.11.24, Big MAP kinase 1, BMK-1, Extracellular signal-regulated kinase 5, ERK-5, Mapk7, Bmk1, Erk5
Accession No.	Q9WVS8
Target Species	Mus musculus (Mouse)
Calculated MW	36.17 kDa
Target Length	Glu56-Asp352
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.