

Recombinant Mouse LOXL4 Protein, N-His

Catalog No: #AP97925

Package Size: #AP97925-1 50ug #AP97925-2 100ug #AP97925-3 1mg

Description

Product Name	Recombinant Mouse LOXL4 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Lysyl oxidase homolog 4, 1.4.3.13, Lysyl oxidase-like protein 4, Lysyl oxidase-related protein C, Loxl4, Loxc
Accession No.	Q924C6
Target Species	Mus musculus (Mouse)
Calculated MW	45.18 kDa
Target Length	Ser143-Asn531
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.