

Recombinant ASFV NP868R Protein, N-His

Catalog No: #AP97950

Package Size: #AP97950-1 50ug #AP97950-2 100ug #AP97950-3 1mg

Description

Product Name	Recombinant ASFV NP868R Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	mRNA-capping enzyme, VTF/CE, Polynucleotide 5'-triphosphatase, 3.6.1.74, mRNA 5'-triphosphatase, TPase, mRNA guanylyltransferase, 2.7.7.50, GTP--RNA guanylyltransferase, GTase, mRNA (guanine-N(7))-methyltransferase, 2.1.1.56, Ba71V-101, NP868R
Accession No.	P32094
Target Species	African swine fever virus (ASFV)
Calculated MW	34.81 kDa
Target Length	Met586-Asn868
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.