

Recombinant ASFV EP402R Protein, N-His

Catalog No: #AP97960

Package Size: #AP97960-1 50ug #AP97960-2 100ug #AP97960-3 1mg

Description

Product Name	Recombinant ASFV EP402R Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	CD2 homolog, 5HL, CD2v, T-lymphocyte CD2 receptor-like protein, pEP402R, EP402R, EP402R CDS, ASFV-Georgia_4-082, ASFV_Kyiv_2016_131_00086
Accession No.	A0A0A7W5A6
Target Species	African swine fever virus (ASFV)
Calculated MW	17.11 kDa
Target Length	Ser230-Ile360
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.