

Recombinant ASFV EP153R Protein, N-His

Catalog No: #AP97961

Package Size: #AP97961-1 50ug #AP97961-2 100ug #AP97961-3 1mg

Description

Product Name	Recombinant ASFV EP153R Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Lectin-like protein EP153R, EP153R, EP153R CDS, ASFV-Georgia_4-081, ASFV_Kyiv_2016_131_00085
Accession No.	A0A0A7W5A0
Target Species	African swine fever virus (ASFV)
Calculated MW	15.64 kDa
Target Length	Gly46-Lys158
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.