

Recombinant Mouse Ins1/Insulin-1 Protein, C-Fc

Catalog No: #AP97963

Package Size: #AP97963-1 50ug #AP97963-2 100ug #AP97963-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Ins1/Insulin-1 Protein, C-Fc
Host Species	Insect Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Insulin-1, Insulin-1 B chain, Insulin-1 A chain, Ins1, Ins-1
Accession No.	P01325
Target Species	Mus musculus (Mouse)
Calculated MW	38.19 kDa
Target Length	Met1-Asn108
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.