

Recombinant Streptococcus suis ApuA Protein, N-His

Catalog No: #AP97985

Package Size: #AP97985-1 50ug #AP97985-2 100ug #AP97985-3 1mg

Description

Product Name	Recombinant Streptococcus suis ApuA Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Pullulanase, 3.2.1.41, Alpha-dextrin endo-1,6-alpha-glucosidase, Pullulan 6-glucanohydrolase, SSU05_2064
Accession No.	A4VY41
Target Species	Streptococcus suis
Calculated MW	116.56 kDa
Target Length	Ile57-Thr1098
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.