

Recombinant Common carp CYP19A1 Protein, C-His

Catalog No: #AP97998

Package Size: #AP97998-1 50ug #AP97998-2 100ug #AP97998-3 1mg

Description

Product Name	Recombinant Common carp CYP19A1 Protein, C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	ovarian aromatase isoform X2, LOC109072133
Accession No.	XP_042630895.1
Target Species	Common carp
Calculated MW	51.20 kDa
Target Length	Ser88-Leu518
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.