

Recombinant *Neisseria meningitidis* NadA Protein, N-His

Catalog No: #AP98004

Package Size: #AP98004-1 50ug #AP98004-2 100ug #AP98004-3 1mg

Description

Product Name	Recombinant <i>Neisseria meningitidis</i> NadA Protein, N-His
Host Species	<i>E. coli</i>
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	autotransporter adhesin NadA, <i>Neisseria</i> adhesin A, NadA
Accession No.	WP_186003383.1
Target Species	<i>Neisseria meningitidis</i>
Calculated MW	16.36 kDa
Target Length	Ala39-Ala170
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.