

## Recombinant Human EP300 Protein, N-His

Catalog No: #AP98016

Package Size: #AP98016-1 50ug #AP98016-2 100ug #AP98016-3 1mg

## Description

|                |                                                                                                                                                                                                                                                                     |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant Human EP300 Protein, N-His                                                                                                                                                                                                                              |
| Host Species   | E. coli                                                                                                                                                                                                                                                             |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                                                                                                                     |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                                                                                                                     |
| Other Names    | Histone crotonyltransferase p300, Protein lactyltransferase p300, Protein propionyltransferase p300, Histone acetyltransferase p300, Protein 2-hydroxyisobutyltransferase p300, EP300, p300 HAT, Histone butyryltransferase p300, P300, E1A-associated protein p300 |
| Accession No.  | Q09472                                                                                                                                                                                                                                                              |
| Target Species | Homo sapiens (Human)                                                                                                                                                                                                                                                |
| Calculated MW  | 31.76 kDa                                                                                                                                                                                                                                                           |
| Target Length  | Glu1351-Ile1601                                                                                                                                                                                                                                                     |
| Tag Info       | His-tag                                                                                                                                                                                                                                                             |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                                                                                                          |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.                                                                                           |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.