

Recombinant HAdV-2 E1B-19K Protein, N-GST & C-His

Catalog No: #AP98019

Package Size: #AP98019-1 50ug #AP98019-2 100ug #AP98019-3 1mg

Description

Product Name	Recombinant HAdV-2 E1B-19K Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	E1B protein, small T-antigen, E1B 19 kDa protein, E1B-19K, E1B-175R, E1B
Accession No.	P03247
Target Species	Human adenovirus C serotype 2 (HAdV-2) (Human adenovirus 2)
Calculated MW	36.67 kDa
Target Length	Ser35-His108
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.