

## Recombinant Human PCMT1 Protein, N-His

Catalog No: #AP98035

Package Size: #AP98035-1 50ug #AP98035-2 100ug #AP98035-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                |                                                                                                                                                                                                                            |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant Human PCMT1 Protein, N-His                                                                                                                                                                                     |
| Host Species   | E. coli                                                                                                                                                                                                                    |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                                                                            |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                                                                            |
| Other Names    | PIMT, L-isoaspartyl protein carboxyl methyltransferase, PCMT1, Protein-L-isoaspartate(D-aspartate)<br>O-methyltransferase, Protein L-isoaspartyl/D-aspartyl methyltransferase, Protein-beta-aspartate<br>methyltransferase |
| Accession No.  | P22061                                                                                                                                                                                                                     |
| Target Species | Homo sapiens (Human)                                                                                                                                                                                                       |
| Calculated MW  | 26.80 kDa                                                                                                                                                                                                                  |
| Target Length  | Met1-Lys227                                                                                                                                                                                                                |
| Tag Info       | His-tag                                                                                                                                                                                                                    |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                                                                 |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store<br>at -20 to -80°C for twelve months from the date of receipt.                                               |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.