

Recombinant Human POMC Protein, N-GST & C-His

Catalog No: #AP98037

Package Size: #AP98037-1 50ug #AP98037-2 100ug #AP98037-3 1mg

Description

Product Name	Recombinant Human POMC Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Melanotropin alpha, Beta-MSH, Gamma-MSH, Beta-LPH, Melanotropin beta, Adrenocorticotrophic hormone, Pro-opiomelanocortin, ACTH, CLIP, Gamma-LPH, Corticotropin-lipotropin, Alpha-MSH, POMC
Accession No.	P01189
Target Species	Homo sapiens (Human)
Calculated MW	39.18 kDa
Target Length	Ala169-Glu267
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.