

Recombinant EBV/HHV-4 BZLF1 Protein, N-GST & C-His

Catalog No: #AP98038

Package Size: #AP98038-1 50ug #AP98038-2 100ug #AP98038-3 1mg

Description

Product Name	Recombinant EBV/HHV-4 BZLF1 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Lytic switch protein BZLF1, EB1, Protein Z, Trans-activator protein BZLF1, Zta, bZIP transcription factor ZEBRA, BZLF1
Accession No.	P03206
Target Species	Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)
Calculated MW	36.56 kDa
Target Length	Leu175-Phe245
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.