

Recombinant Mouse PDHA1 Protein, N-His

Catalog No: #AP98040

Package Size: #AP98040-1 50ug #AP98040-2 100ug #AP98040-3 1mg

Description

Product Name	Recombinant Mouse PDHA1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial, 1.2.4.1, PDHE1-A type I, Pdha1, Pdha-1
Accession No.	P35486
Target Species	Mus musculus (Mouse)
Calculated MW	43 kDa
Target Length	Phe30-Ser390
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.