

Recombinant ZIKV Envelope protein E1 Protein, N-His

Catalog No: #AP98062

Package Size: #AP98062-1 50ug #AP98062-2 100ug #AP98062-3 1mg

Description

Product Name	Recombinant ZIKV Envelope protein E1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Genome polyprotein, Envelope protein E
Accession No.	A0A024B7W1
Target Species	Zika virus (isolate ZIKV/Human/French Polynesia/10087PF/2013) (ZIKV)
Calculated MW	52.14 kDa
Target Length	Ile291-Ser751
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.