

Recombinant Human PDIA4 Protein, C-His

Catalog No: #AP98065

Package Size: #AP98065-1 50ug #AP98065-2 100ug #AP98065-3 1mg

Description

Product Name	Recombinant Human PDIA4 Protein, C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	ERp-72, ER protein 70, ER protein 72, ERP70, Protein disulfide-isomerase A4, ERP72, PDIA4, Endoplasmic reticulum resident protein 72, ERp72, ERp70, Endoplasmic reticulum resident protein 70
Accession No.	P13667
Target Species	Homo sapiens (Human)
Calculated MW	71.86 kDa
Target Length	Val21-Leu645
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.