

Recombinant Human PTPRN Protein, N-MAT

Catalog No: #AP98069

Package Size: #AP98069-1 50ug #AP98069-2 100ug #AP98069-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PTPRN Protein, N-MAT
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	ICA512-TMF, ICA512-CCF, Receptor-type tyrosine-protein phosphatase-like N, Islet cell autoantigen 3, ICA512, PTPRN, Islet cell antigen 512, ICA 512, PTP IA-2, ICA3, ICA512-NTF, R-PTP-N
Accession No.	Q16849
Target Species	Homo sapiens (Human)
Calculated MW	44.79 kDa
Target Length	His603-Gln979
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.