

Recombinant Human MAGEA6 Protein, N-His

Catalog No: #AP98073

Package Size: #AP98073-1 50ug #AP98073-2 100ug #AP98073-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MAGEA6 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	MAGEA6, Cancer/testis antigen 1.6, MAGE6, MAGE3B antigen, MAGE-6 antigen, CT1.6, Melanoma-associated antigen 6
Accession No.	P43360
Target Species	Homo sapiens (Human)
Calculated MW	27.19 kDa
Target Length	Ser96-Glu314
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.