

Recombinant Human MCCC2 Protein, N-His

Catalog No: #AP98074

Package Size: #AP98074-1 50ug #AP98074-2 100ug #AP98074-3 1mg

Description

Product Name	Recombinant Human MCCC2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	3-methylcrotonyl-CoA carboxylase non-biotin-containing subunit, 3-methylcrotonyl-CoA:carbon dioxide ligase subunit beta, Methylcrotonoyl-CoA carboxylase beta chain, mitochondrial, MCCB, MCCC2, 3-methylcrotonyl-CoA carboxylase 2, MCCase subunit beta
Accession No.	Q9HCC0
Target Species	Homo sapiens (Human)
Calculated MW	25.67 kDa
Target Length	Lys80-Asn295
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.