

Recombinant Human NEK2 Protein, N-His

Catalog No: #AP98077

Package Size: #AP98077-1 50ug #AP98077-2 100ug #AP98077-3 1mg

Description

Product Name	Recombinant Human NEK2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	NLK1, NEK2, NimA-like protein kinase 1, HSPK 21, NEK2A, Serine/threonine-protein kinase Nek2, NimA-related protein kinase 2, Never in mitosis A-related kinase 2
Accession No.	P51955
Target Species	Homo sapiens (Human)
Calculated MW	33.71 kDa
Target Length	Met1-Ile271
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.