

Recombinant EBV/HHV-4 gH/BXLF2 Protein, C-His

Catalog No: #AP98080

Package Size: #AP98080-1 50ug #AP98080-2 100ug #AP98080-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant EBV/HHV-4 gH/BXLF2 Protein, C-His
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Envelope glycoprotein H, gH, gH, BXLF2
Accession No.	P03231
Target Species	Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)
Calculated MW	37.66 kDa
Target Length	Met361-Tyr674
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.