

Recombinant CHIKV Spike glycoprotein E2 Protein, C-His

Catalog No: #AP98083

Package Size: #AP98083-1 50ug #AP98083-2 100ug #AP98083-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant CHIKV Spike glycoprotein E2 Protein, C-His
Host Species	Mammalian
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	structural polyprotein, structural polyprotein
Accession No.	AIL00833.1
Target Species	Chikungunya virus (CHIKV)
Calculated MW	41.26 kDa
Target Length	Ser326-Gln666
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.