

Recombinant Human ZP3 Protein, N-His

Catalog No: #AP98090

Package Size: #AP98090-1 50ug #AP98090-2 100ug #AP98090-3 1mg

Description

Product Name	Recombinant Human ZP3 Protein, N-His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	ZP3B, Zona pellucida protein C, ZP3, Zp-3, ZP3A/ZP3B, Sperm receptor, ZPC, Zona pellucida glycoprotein 3, ZP3A, Zona pellucida sperm-binding protein 3
Accession No.	P21754
Target Species	Homo sapiens (Human)
Calculated MW	43.82 kDa
Target Length	Gln23-Val387
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.