

Recombinant Mouse TENT4B Protein, N-His

Catalog No: #AP98091

Package Size: #AP98091-1 50ug #AP98091-2 100ug #AP98091-3 1mg

Description

Product Name	Recombinant Mouse TENT4B Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Terminal nucleotidyltransferase 4B, Non-canonical poly(A) RNA polymerase PAPD5, 2.7.7.19, PAP-associated domain-containing protein 5, Terminal guanylyltransferase, 2.7.7.-, Terminal uridylyltransferase 3, TUTase 3, Topoisomerase-related function protein 4-2, TRF4-2, Tent4b, Papd5
Accession No.	Q68ED3
Target Species	Mus musculus (Mouse)
Calculated MW	70.81 kDa
Target Length	Phe2-Arg633
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.